

SRP9 Polyclonal Antibody

Catalog No: #30522



Package Size: #30522-1 50ul #30522-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SRP9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human SRP9 (NP_003124.1).
Other Names	SRP9; ALURBP; signal recognition particle 9
Accession No.	Swiss-Prot#:P49458NCBI Gene ID:6726
Uniprot	P49458
GeneID	6726;
Calculated MW	12kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:1000 - 1:2000IF□1:50 - 1:200

Images



Western blot analysis of extracts of A-431 cells, using SRP9 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only