

PIGQ Antibody

Catalog No: #46128



Package Size: #46128-1 50ul #46128-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PIGQ Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PIGQ Antibody detects endogenous levels of total PIGQ
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PIGQ
Target Name	PIGQ
Other Names	PIGQ, C407A10.1, PIG-Q, GPI1, HGPI1
Accession No.	Swiss-Prot#: Q9BRB3NCBI Gene ID:9091
Uniprot	Q9BRB3
GeneID	9091;
Calculated MW	84 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only