

PLAA Polyclonal Antibody

Catalog No: #30542



Package Size: #30542-1 50ul #30542-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PLAA Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human PLAA (NP_001026859.1).
Other Names	PLAA; DOA1; PLA2P; PLAP; NDMSBA; phospholipase A2 activating protein
Accession No.	Swiss-Prot#:Q9Y263NCBI Gene ID:9373
Uniprot	Q9Y263
GeneID	9373;
Calculated MW	87kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PLAA at 1:1000 dilution.

Background

Note: This product is for in vitro research use only