

GOSR1 Polyclonal Antibody

Catalog No: #30545

Package Size: #30545-1 50ul #30545-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GOSR1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human GOSR1 (NP_004862.1).
Other Names	GOSR1;GOLIM2;GOS-28;GOS28;GOS28/P28;GS28;P28
Accession No.	Uniprot:O95249GeneID:9527
Uniprot	O95249
GeneID	9527
Calculated MW	29kDa
SDS-PAGE MW	29kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB □ 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse testis, using GOSR1 antibody.

Background

This gene encodes a trafficking membrane protein which transports proteins among the endoplasmic reticulum and the Golgi and between Golgi compartments. This protein is considered an essential component of the Golgi SNAP receptor (SNARE) complex. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

Note: This product is for in vitro research use only