

RABEPK Polyclonal Antibody

Catalog No: #30557



Package Size: #30557-1 50ul #30557-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RABEPK Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human RABEPK (NP_005824.2).
Other Names	RABEPK; RAB9P40; bA65N13.1; p40; rab9 effector protein with kelch motifs
Accession No.	Swiss-Prot#:Q7Z6M1NCBI Gene ID:10244
Uniprot	Q7Z6M1
GeneID	10244;
Calculated MW	45kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of HeLa cells, using RABEPK at 1:3000 dilution.

Background

Note: This product is for in vitro research use only