

PNMA2 Polyclonal Antibody

Catalog No: #30567



Package Size: #30567-1 50ul #30567-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PNMA2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Rat
Immunogen Description	Recombinant fusion protein of human PNMA2 (NP_009188.1).
Other Names	PNMA2; MA2; MM2; RGAG2; paraneoplastic antigen Ma2
Accession No.	Swiss-Prot#:Q9UL42NCBI Gene ID:10687
Uniprot	Q9UL42
GeneID	10687;
Calculated MW	41kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PNMA2 at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021). Exposure time: 90s.

Background

Note: This product is for in vitro research use only