

LARP1 Polyclonal Antibody

Catalog No: #30580



Package Size: #30580-1 50ul #30580-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LARP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant Protein of human LARP1
Other Names	LARP1; LARP; la-related protein 1
Accession No.	Swiss-Prot#:Q6PKG0NCBI Gene ID:23367
Uniprot	Q6PKG0
GeneID	23367;
Calculated MW	163kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using LARP1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only