

RBM15 Polyclonal Antibody

Catalog No: #30623



Package Size: #30623-1 50ul #30623-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RBM15 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human RBM15 (NP_001188474.1).
Other Names	RBM15; OTT; OTT1; SPEN; RNA binding motif protein 15
Accession No.	Swiss-Prot#:Q96T37NCBI Gene ID:64783
Uniprot	Q96T37
GeneID	64783;
Calculated MW	107kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using RBM15 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only