

LPCAT1 Polyclonal Antibody

Catalog No: #30627



Package Size: #30627-1 50ul #30627-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LPCAT1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human LPCAT1 (NP_079106.3).
Other Names	LPCAT1; AGPAT10; AGPAT9; AYTL2; LPCAT-1; PFAAP3; lpcat; lysoPAFAT; lysophosphatidylcholine acyltransferase 1
Accession No.	Swiss-Prot#:Q8NF37NCBI Gene ID:79888
Uniprot	Q8NF37
GeneID	79888;
Calculated MW	59kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:1000 - 1:2000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using LPCAT1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only