

RMI1 Polyclonal Antibody

Catalog No: #30629

Package Size: #30629-1 50ul #30629-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RMI1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human RMI1 (NP_079221.2).
Other Names	RMI1; BLAP75; C9orf76; FAAP75; RecQ mediated genome instability 1
Accession No.	Swiss-Prot#:Q9H9A7NCBI Gene ID:80010
Uniprot	Q9H9A7
GeneID	80010;
Calculated MW	75kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using RMI1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only