

GRPEL1 Polyclonal Antibody

Catalog No: #30632



Package Size: #30632-1 50ul #30632-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GRPEL1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human GRPEL1 (NP_079472.1).
Other Names	GRPEL1; GrpE; HMGE; mt-GrpE#1; GrpE like 1, mitochondrial
Accession No.	Swiss-Prot#:Q9HAV7NCBI Gene ID:80273
Uniprot	Q9HAV7
GeneID	80273;
Calculated MW	20kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using GRPEL1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only