

MTDH Polyclonal Antibody

Catalog No: #30670

Package Size: #30670-1 50ul #30670-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MTDH Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human MTDH (NP_848927.2).
Other Names	MTDH; 3D3; AEG-1; AEG1; LYRIC; LYRIC/3D3; metadherin
Accession No.	Swiss-Prot#:Q86UE4NCBI Gene ID:92140
Uniprot	Q86UE4
GeneID	92140;
Calculated MW	64kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200IF□1:50 - 1:200


Images




Immunohistochemistry of paraffin-embedded human gastric using MTDH at dilution of 1:200 (40x lens).




Immunohistochemistry of paraffin-embedded mouse brain using MTDH at dilution of 1:200 (40x lens).



Immunofluorescence analysis of U2OS cells using MTDH at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunohistochemistry of paraffin-embedded rat brain using MTDH at dilution of 1:200 (40x lens).



Western blot analysis of extracts of 22Rv1 cells, using MTDH at 1:1000 dilution.

Background

Note: This product is for in vitro research use only