

PWP2 Polyclonal Antibody

Catalog No: #30678



Package Size: #30678-1 50ul #30678-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PWP2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human PWP2 (NP_005040.2).
Other Names	PWP2; EHOC-17; PWP2H; UTP1; periodic tryptophan protein 2 homolog
Accession No.	Swiss-Prot#:Q15269NCBI Gene ID:5822
Uniprot	Q15269
GeneID	5822;
Calculated MW	102kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:200 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using PWP2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only