

ZCCHC11 Polyclonal Antibody

Catalog No: #30686



Package Size: #30686-1 50ul #30686-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZCCHC11 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZCCHC11 (NP_001009881.1).
Other Names	ZCCHC11;PAPD3;TUT4
Accession No.	Uniprot:Q5TAX3GeneID:23318
Uniprot	Q5TAX3
GeneID	23318
Calculated MW	185kDa
SDS-PAGE MW	185kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:1000IF 1:50 - 1:200


Images




Immunofluorescence analysis of U2OS cells using ZCCHC11 Rabbit pAb.



Western blot analysis of extracts of 293T cells, using ZCCHC11 Antibody.



Immunofluorescence analysis of H9C2 cells using ZCCHC11 Rabbit pAb.



Immunofluorescence analysis of L929 cells using ZCCHC11 Rabbit pAb.

Background

ZCCHC11 is an RNA uridylyltransferase (EC 2.7.7.52) that uses UTP to add uridines to the 3-prime end of substrate RNA molecules (Jones et al., 2009 [PubMed 19701194]).

Note: This product is for in vitro research use only