

TDRD3 Polyclonal Antibody

Catalog No: #30690



Package Size: #30690-1 50ul #30690-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TDRD3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human TDRD3 (NP_110421.1).
Other Names	TDRD3; tudor domain containing 3
Accession No.	Swiss-Prot#:Q9H7E2NCBI Gene ID:81550
Uniprot	Q9H7E2
GeneID	81550;
Calculated MW	73kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of U2OS cells, using TDRD3 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only