

AQR Polyclonal Antibody

Catalog No: #30691

Package Size: #30691-1 50ul #30691-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	AQR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human AQR (NP_055506.1).
Other Names	AQR; IBP160; fSAP164; intron-binding protein aquarius
Accession No.	Swiss-Prot#:O60306NCBI Gene ID:9716
Uniprot	O60306
GeneID	9716;
Calculated MW	171kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB□1:500 - 1:1000IHC□1:50 - 1:100


Images




Immunohistochemistry of paraffin-embedded human lung cancer using AQR at dilution of 1:100 (40x lens).



Western blot analysis of extracts of various cell lines, using AQR at 1:3000 dilution.



Immunohistochemistry of paraffin-embedded mouse testis using AQR at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded rat spleen using AQR at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only