

RBBP6 Polyclonal Antibody

Catalog No: #30775

Package Size: #30775-1 50ul #30775-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RBBP6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human RBBP6 (NP_116015.2).
Other Names	RBBP6; MY038; P2P-R; PACT; RBQ-1; SNAMA; E3 ubiquitin-protein ligase RBBP6
Accession No.	Swiss-Prot#:Q7Z6E9NCBI Gene ID:5930
Uniprot	Q7Z6E9
GeneID	5930;
Calculated MW	160kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

IF□1:50 - 1:200

Images



Immunofluorescence analysis of A-549 cells using RBBP6 .
Blue: DAPI for nuclear staining.

Background

The retinoblastoma tumor suppressor (pRB) protein binds with many other proteins. In various human cancers, pRB suppresses cellular proliferation and is inactivated. Cell cycle-dependent phosphorylation regulates the activity of pRB. This gene encodes a protein which binds to underphosphorylated but not phosphorylated pRB. Multiple alternatively spliced transcript variants that encode different isoforms have been found for this gene.

Note: This product is for in vitro research use only