

## PPP1R3D Antibody

Catalog No: #46153

Package Size: #46153-1 50ul #46153-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	PPP1R3D Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PPP1R3D Antibody detects endogenous levels of total PPP1R3D
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PPP1R3D
Target Name	PPP1R3D
Other Names	PPP1R6
Accession No.	Swiss-Prot#: O95685NCBI Gene ID:5509
Uniprot	O95685
GeneID	5509;
Calculated MW	33 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4o½ over night

---

Note: This product is for in vitro research use only