

PPM1E Antibody

Catalog No: #46154



Package Size: #46154-1 50ul #46154-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PPM1E Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	PPM1E Antibody detects endogenous levels of total PPM1E
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PPM1E
Target Name	PPM1E
Other Names	CAMKN, KIAA1072, POPX1
Accession No.	Swiss-Prot#: Q8WY54NCBI Gene ID:22843
Uniprot	Q8WY54
GeneID	22843;
Calculated MW	85 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of Hela lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only