

KIAA1456 Polyclonal Antibody

Catalog No: #30845



Package Size: #30845-1 50ul #30845-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	KIAA1456 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human KIAA1456 (NP_065895.2).
Other Names	KIAA1456;C8orf79;TRM9L
Accession No.	Uniprot:Q9P272GeneID:57604
Uniprot	Q9P272
GeneID	57604
Calculated MW	51kDa
SDS-PAGE MW	51kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using KIAA1456 antibody.

Background

May modifie wobble uridines in specific arginine and glutamic acid tRNAs. Acts as a tumor suppressor by promoting the expression of LIN9.

Note: This product is for in vitro research use only