

KLF3 Polyclonal Antibody

Catalog No: #30846



Package Size: #30846-1 50ul #30846-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	KLF3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human KLF3 (NP_057615.3).
Other Names	KLF3; BKLF; Kruppel like factor 3
Accession No.	Swiss-Prot#:P57682NCBI Gene ID:51274
Uniprot	P57682
GeneID	51274;
Calculated MW	39kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100

Images



Western blot analysis of extracts of various cell lines, using KLF3 at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021). Exposure time: 60s.



Immunofluorescence analysis of A549 cells using KLF3 .

Background

Note: This product is for in vitro research use only