

NAT8B Polyclonal Antibody

Catalog No: #30852



Package Size: #30852-1 50ul #30852-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NAT8B Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NAT8B (NP_057431.2).
Other Names	NAT8B; CML2; Hcml2; NAT8BP; putative N-acetyltransferase 8B
Accession No.	Swiss-Prot#:Q9UHF3NCBI Gene ID:51471
Uniprot	Q9UHF3
GeneID	51471;
Calculated MW	25kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200


Images




Western blot analysis of extracts of rat brain, using NAT8B at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded rat kidney using NAT8B at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded human stomach using NAT8B at dilution of 1:200 (40x lens).



Immunohistochemistry of paraffin-embedded mouse heart using NAT8B at dilution of 1:200 (40x lens).

Background

Note: This product is for in vitro research use only