

PECR Polyclonal Antibody

Catalog No: #30855



Package Size: #30855-1 50ul #30855-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PECR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Mouse
Immunogen Description	Recombinant fusion protein of human PECR (NP_060911.2).
Other Names	PECR; DCCRPP; HPDHASE; HSA250303; PVIARL; SDR29C1; TERP; peroxisomal trans-2-enoyl-CoA reductase
Accession No.	Swiss-Prot#:Q9BY49NCBI Gene ID:55825
Uniprot	Q9BY49
GeneID	55825;
Calculated MW	33kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200

Images



Western blot analysis of extracts of mouse kidney, using PECR at 1:1000 dilution.

Background

Note: This product is for in vitro research use only