

PHIP Polyclonal Antibody

Catalog No: #30856



Package Size: #30856-1 50ul #30856-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PHIP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human PHIP (NP_060404.3).
Other Names	PHIP; BRWD2; DCAF14; WDR11; ndrp; PH-interacting protein
Accession No.	Swiss-Prot#:Q8WWQ0NCBI Gene ID:55023
Uniprot	Q8WWQ0
GeneID	55023;
Calculated MW	250kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:100


Images




Immunohistochemistry of paraffin-embedded human well-differentiated squamous skin carcinoma using PHIP at dilution of 1:100 (40x lens).




Immunohistochemistry of paraffin-embedded human kidney cancer using PHIP at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse tumor using PHIP at dilution of 1:100 (40x lens).



Immunofluorescence analysis of MCF-7 cells using PHIP .



Western blot analysis of extracts of various cell lines, using PHIP at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021). Exposure time: 90s.

Background

Note: This product is for in vitro research use only