

PHYHD1 Polyclonal Antibody

Catalog No: #30857



Package Size: #30857-1 50ul #30857-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PHYHD1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PHYHD1 (NP_001094346.1).
Other Names	PHYHD1; phytanoyl-CoA dioxygenase domain containing 1
Accession No.	Swiss-Prot#:Q5SRE7NCBI Gene ID:254295
Uniprot	Q5SRE7
GeneID	254295;
Calculated MW	32kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100

Images



Immunofluorescence analysis of A-549 cells using PHYHD1 .
Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using PHYHD1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only