

ZDHHC16 Polyclonal Antibody

Catalog No: #30865



Package Size: #30865-1 50ul #30865-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZDHHC16 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZDHHC16 (NP_932162.1).
Other Names	ZDHHC16; APH2; DHHC-16; zinc finger DHHC-type containing 16
Accession No.	Swiss-Prot#:Q969W1NCBI Gene ID:84287
Uniprot	Q969W1
GeneID	84287;
Calculated MW	44kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZDHHC16 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only