

C10orf32 Polyclonal Antibody

Catalog No: #30877

Package Size: #30877-1 50ul #30877-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	C10orf32 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human C10orf32 (NP_653192.2).
Other Names	BORCS7; C10orf32; BLOC-1-related complex subunit 7
Accession No.	Swiss-Prot#:Q96B45NCBI Gene ID:119032
Uniprot	Q96B45
GeneID	119032;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:100

Images



Immunofluorescence analysis of MCF-7 cells using C10orf32 .
Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only