

SETD5 Polyclonal Antibody

Catalog No: #30884

Package Size: #30884-1 50ul #30884-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SETD5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human SETD5 (NP_001073986.1).
Other Names	SETD5; SET domain containing 5
Accession No.	Swiss-Prot#:Q9C0A6NCBI Gene ID:55209
Uniprot	Q9C0A6
GeneID	55209;
Calculated MW	158kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200IF□1:50 - 1:100

Images



Western blot analysis of extracts of MCF-7 cells, using SETD5 at 1:300 dilution.

Background

This function of this gene has yet to be determined but mutations in this gene have been associated with autosomal dominant mental retardation-23. Alternative splicing results in multiple transcript variants encoding different isoforms.

Note: This product is for in vitro research use only