

TRIAP1 Polyclonal Antibody

Catalog No: #30887



Package Size: #30887-1 50ul #30887-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TRIAP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human TRIAP1 (NP_057483.1).
Other Names	TRIAP1; HSPC132; MDM35; P53CSV; WF-1; TP53-regulated inhibitor of apoptosis 1
Accession No.	Swiss-Prot#:O43715NCBI Gene ID:51499
Uniprot	O43715
GeneID	51499;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100

Images



Immunofluorescence analysis of MCF-7 cells using TRIAP1 .
Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only