

YTHDC1 Polyclonal Antibody

Catalog No: #30888



Package Size: #30888-1 50ul #30888-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	YTHDC1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human YTHDC1 (NP_001026902.1).
Other Names	YTHDC1; YT521; YT521-B; YTH domain containing 1
Accession No.	Swiss-Prot#:Q96MU7NCBI Gene ID:91746
Uniprot	Q96MU7
GeneID	91746;
Calculated MW	118kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100

Images



Immunofluorescence analysis of A-549 cells using YTHDC1 .
Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using YTHDC1 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only