

ZNF20 Polyclonal Antibody

Catalog No: #30895



Package Size: #30895-1 50ul #30895-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF20 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZNF20 (NP_066966.2).
Other Names	ZNF20; KOX13; zinc finger protein 20
Accession No.	Swiss-Prot#:P17024NCBI Gene ID:7568
Uniprot	P17024
GeneID	7568;
Calculated MW	62kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using ZNF20 at 1:1000 dilution.



Immunohistochemistry of paraffin-embedded human brain cancer using ZNF20 at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded human kidney cancer using ZNF20 at dilution of 1:100 (40x lens).

Background

Note: This product is for in vitro research use only