

DPF1 Polyclonal Antibody

Catalog No: #30897

Package Size: #30897-1 50ul #30897-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DPF1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human DPF1 (NP_001128627.1).
Other Names	DPF1; BAF45b; NEUD4; neuro-d4; double PHD fingers 1
Accession No.	Swiss-Prot#:Q92782NCBI Gene ID:8193
Uniprot	Q92782
GeneID	8193;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000IF□1:50 - 1:100

Images



Immunofluorescence analysis of A549 cells using DPF1 .
Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only