

ZNF264 Polyclonal Antibody

Catalog No: #30901

Package Size: #30901-1 50ul #30901-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ZNF264 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Rat
Immunogen Description	Recombinant fusion protein of human ZNF264 (NP_003408.1).
Other Names	ZNF264; zinc finger protein 264
Accession No.	Swiss-Prot#:O43296NCBI Gene ID:9422
Uniprot	O43296
GeneID	9422;
Calculated MW	71kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of rat serum, using ZNF264 at 1:1000 dilution.

Background

This gene encodes a zinc finger protein and belongs to the krueppel C2H2-type zinc-finger protein family. Zinc finger proteins are often localized in the nucleus, bind nucleic acids, and regulate transcription.

Note: This product is for in vitro research use only