

ARHGAP44 Polyclonal Antibody

Catalog No: #30902



Package Size: #30902-1 50ul #30902-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ARHGAP44 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ARHGAP44 (NP_055674.4).
Other Names	ARHGAP44; NPC-A-10; RICH2; rho GTPase-activating protein 44
Accession No.	Swiss-Prot#:Q17R89NCBI Gene ID:9912
Uniprot	Q17R89
GeneID	9912;
Calculated MW	110kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ARHGAP44 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only