

NAA40 Polyclonal Antibody

Catalog No: #30910



Package Size: #30910-1 50ul #30910-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	NAA40 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NAA40 (NP_079047.2).
Other Names	NAA40; NAT11; NatD; PATT1; hNatD; N-alpha-acetyltransferase 40
Accession No.	Swiss-Prot#:Q86UY6NCBI Gene ID:79829
Uniprot	Q86UY6
GeneID	79829;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IF 1:50 - 1:100

Images



Immunofluorescence analysis of A549 cells using NAA40 .
Blue: DAPI for nuclear staining.

Background

Note: This product is for in vitro research use only