

SPIB Polyclonal Antibody

Catalog No: #30928



Package Size: #30928-1 50ul #30928-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SPIB Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human
Immunogen Description	Recombinant fusion protein of human SPIB (NP_003112.2).
Other Names	SPIB;SPI-B
Accession No.	Uniprot:Q01892GeneID:6689
Uniprot	Q01892
GeneID	6689
Calculated MW	46kDa
SDS-PAGE MW	46kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of Raji cells, using SPIB antibody.

Background

The protein encoded by this gene is a transcriptional activator that binds to the PU-box (5'-GAGGAA-3') and acts as a lymphoid-specific enhancer. Four transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only