

ZNF169 Polyclonal Antibody

Catalog No: #30945



Package Size: #30945-1 50ul #30945-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZNF169 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ZNF169 (NP_919301.2).
Other Names	ZNF169
Accession No.	Uniprot:Q14929GeneID:169841
Uniprot	Q14929
GeneID	169841
Calculated MW	72kDa
SDS-PAGE MW	72kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZNF169 antibody.

Background

May be involved in transcriptional regulation.

Note: This product is for in vitro research use only