

ELF2 Polyclonal Antibody

Catalog No: #30946



Package Size: #30946-1 50ul #30946-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ELF2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ELF2 (NP_006865.1).
Other Names	ELF2; EU32; NERF; NERF-1A; NERF-1B; NERF-1a; b; NERF-2; E74 like ETS transcription factor 2
Accession No.	Swiss-Prot#:Q15723NCBI Gene ID:1998
Uniprot	Q15723
GeneID	1998;
Calculated MW	64kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ELF2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only