

ZBTB25 Polyclonal Antibody

Catalog No: #30951



Package Size: #30951-1 50ul #30951-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ZBTB25 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human ZBTB25 (NP_008908.2).
Other Names	ZBTB25;C14orf51;KUP;ZNF46
Accession No.	Uniprot:P24278GenelD:7597
Uniprot	P24278
GenelD	7597
Calculated MW	49kDa
SDS-PAGE MW	49kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ZBTB25 antibody.

Background

May be involved in transcriptional regulation.

Note: This product is for in vitro research use only