

CRIP1 Polyclonal Antibody

Catalog No: #30960



Package Size: #30960-1 50ul #30960-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CRIP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human CRIP1 (NP_001302.1).
Other Names	CRIP1; CRHP; CRIP; CRP-1; CRP1; cysteine-rich protein 1
Accession No.	Swiss-Prot#:P50238NCBI Gene ID:1396
Uniprot	P50238
GeneID	1396;
Calculated MW	13kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse lung, using CRIP1 at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021). Exposure time: 5s.

Background

Note: This product is for in vitro research use only