

UGT2B10 Polyclonal Antibody

Catalog No: #30967



Package Size: #30967-1 50ul #30967-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	UGT2B10 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human UGT2B10 (NP_001066.1).
Other Names	UGT2B10;UDPGT2B10
Accession No.	Uniprot:P36537GenelD:7365
Uniprot	P36537
GenelD	7365
Calculated MW	65kDa/80kDa
SDS-PAGE MW	65kDa/80kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using UGT2B10 antibody.

Background

UDPGT is of major importance in the conjugation and subsequent elimination of potentially toxic xenobiotics and endogenous compounds.

Note: This product is for in vitro research use only