

PLEKHF2 Polyclonal Antibody

Catalog No: #30980



Package Size: #30980-1 50ul #30980-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PLEKHF2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PLEKHF2 (NP_078889.1).
Other Names	PLEKHF2; EAPF; PHAFIN2; ZFYVE18; pleckstrin homology and FYVE domain containing 2
Accession No.	Swiss-Prot#:Q9H8W4NCBI Gene ID:79666
Uniprot	Q9H8W4
GeneID	79666;
Calculated MW	Refer to figures
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB□1:500 - 1:2000

Background

Note: This product is for in vitro research use only