

UGT3A2 Polyclonal Antibody

Catalog No: #30983



Package Size: #30983-1 50ul #30983-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	UGT3A2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human UGT3A2 (NP_777574.2).
Other Names	UGT3A2; UDP-glucuronosyltransferase 3A2
Accession No.	Swiss-Prot#:Q3SY77NCBI Gene ID:167127
Uniprot	Q3SY77
GeneID	167127;
Calculated MW	75kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using UGT3A2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only