

PLCXD2 Polyclonal Antibody

Catalog No: #30986



Package Size: #30986-1 50ul #30986-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PLCXD2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PLCXD2 (NP_001172035.1).
Other Names	PLCXD2; PI-PLC X domain-containing protein 2
Accession No.	Swiss-Prot#:Q0VAA5NCBI Gene ID:257068
Uniprot	Q0VAA5
GeneID	257068;
Calculated MW	30kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using PLCXD2 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only