

YWHAZ Polyclonal Antibody

Catalog No: #30992

Package Size: #30992-1 50ul #30992-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	YWHAZ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human YWHAZ (NP_001129171.1).
Other Names	YWHAZ; 14-3-3-zeta; HEL-S-3; HEL-S-93; HEL4; KCIP-1; YWHAD; 14-3-3 protein zeta/delta
Accession No.	Swiss-Prot#:P63104NCBI Gene ID:7534
Uniprot	P63104
GeneID	7534;
Calculated MW	28kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C


Application Details

WB□1:500 - 1:2000IHC□1:50 - 1:200IF□1:50 - 1:200


Images



Immunohistochemistry of paraffin-embedded rat brain using YWHAZ at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded mouse brain using YWHAZ at dilution of 1:100 (40x lens).



Western blot analysis of extracts of various cell lines, using YWHAZ at 1:1000 dilution.

Background

This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and sheep orthologs. The encoded protein interacts with IRS1 protein, suggesting a role in regulating insulin sensitivity. Several transcript variants that differ in the 5' UTR but that encode the same protein have been identified for this gene.

Note: This product is for in vitro research use only