

PDE4DIP Polyclonal Antibody

Catalog No: #31322



Package Size: #31322-1 50ul #31322-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PDE4DIP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human PDE4DIP (NP_071754.3).
Other Names	PDE4DIP;CMYA2;MMGL
Accession No.	Uniprot:Q5VU43GeneID:9659
Uniprot	Q5VU43
GeneID	9659
Calculated MW	265kDa
SDS-PAGE MW	265kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of mouse skeletal muscle, using PDE4DIP antibody.

Background

The protein encoded by this gene serves to anchor phosphodiesterase 4D to the Golgi/centrosome region of the cell. Defects in this gene may be a cause of myeloproliferative disorder (MBD) associated with eosinophilia. Several transcript variants encoding different isoforms have been found for this gene.

Note: This product is for in vitro research use only