

ALG1 Polyclonal Antibody

Catalog No: #31348



Package Size: #31348-1 50ul #31348-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ALG1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human ALG1 (NP_061982.3).
Other Names	ALG1;CDG1K;HMAT1;HMT-1;HMT1;MT-1;Mat-1;hMat-1
Accession No.	Uniprot:Q9BT22GeneID:56052
Uniprot	Q9BT22
GeneID	56052
Calculated MW	60kDa
SDS-PAGE MW	60kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using ALG1 antibody.

Background

The enzyme encoded by this gene catalyzes the first mannosylation step in the biosynthesis of lipid-linked oligosaccharides. This gene is mutated in congenital disorder of glycosylation type Ik.

Note: This product is for in vitro research use only