

TM9SF3 Polyclonal Antibody

Catalog No: #31349



Package Size: #31349-1 50ul #31349-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TM9SF3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	Recombinant fusion protein of human TM9SF3 (NP_064508.3).
Other Names	TM9SF3; EP70-P-iso; SMBP; transmembrane 9 superfamily member 3
Accession No.	Swiss-Prot#:Q9HD45NCBI Gene ID:56889
Uniprot	Q9HD45
GeneID	56889;
Calculated MW	90kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using TM9SF3 at 1:1000 dilution.

Background

Note: This product is for in vitro research use only