

REEP1 Polyclonal Antibody

Catalog No: #31353



Package Size: #31353-1 50ul #31353-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	REEP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human REEP1 (NP_075063.1).
Other Names	REEP1;C2orf23;HMN5B;SPG31;Yip2a
Accession No.	Uniprot:Q9H902GeneID:65055
Uniprot	Q9H902
GeneID	65055
Calculated MW	22kDa
SDS-PAGE MW	22kDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using REEP1 antibody.

Background

This gene encodes a mitochondrial protein that functions to enhance the cell surface expression of odorant receptors. Mutations in this gene cause spastic paraplegia autosomal dominant type 31, a neurodegenerative disorder. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only