

NPTN Polyclonal Antibody

Catalog No: #31395

Package Size: #31395-1 50ul #31395-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NPTN Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human NPTN (NP_059429.1).
Other Names	NPTN; GP55; GP65; SDFR1; SDR1; np55; np65; neuroplastin
Accession No.	Swiss-Prot#:Q9Y639NCBI Gene ID:27020
Uniprot	Q9Y639
GeneID	27020;
Calculated MW	45kDa
Formulation	Avoid freeze / thaw cycles. Buffer: PBS with 50% glycerol, pH7.4.
Storage	Store at -20°C

Application Details

WB 1:500 - 1:2000

Images



Western blot analysis of extracts of various cell lines, using NPTN at 1:1000 dilution.

Background

This gene encodes a type I transmembrane protein belonging to the Ig superfamily. The protein is believed to be involved in cell-cell interactions or cell-substrate interactions. Alternative splicing results in multiple transcript variants.

Note: This product is for in vitro research use only